

RESPIRATOR MASK

KN95 tested in accordance with GB2626-2006



ORDER INFORMATION

Order no.	Description	Colour	Sales unit
97365-1000	Respirator mask	white	1.000 pieces (20 x 50 pcs.)

MATERIAL COMPOSITION

- Polypropylene spunbonded fabric
- Elastic bands (2 pcs. around the ears)

DESCRIPTION

- Ergonomically shaped
- The nose clip allows it to be perfectly to the face shape
- High-quality and skin-friendly
- BFE \geq 95%
- Unsterile
- Lint-free
- With two elastic bands around the ears
- Universal size
- Without valve
- Certified according to EN149:2001 + A1:2009

SPECIAL STORAGE INSTRUCTIONS

- First in - first out
- Store in a dry place

USAGE

- Doctors' surgeries
- Hospitals
- Rehab, retirement and care facilities

Expansion approx. (in %):	N/A
Sterilisation:	non-sterile
Washable up to:	N/A
CE mark:	no
Customs tariff no.:	6307.9099
Shelf life:	5 years
Classification:	I
Single-use product:	yes

CERTIFICATION



Last updated: 06-2020